



OEM PARTNER VITAL SIGNS MONITOR CATALOG

Complete list of
VSM, Thermometer,
and Cable options
with part #'s

Table of Contents **Page #**

Introduction	2
Alfamed	3
Capsule	4
Direct Supply	5
Draeger	6
Edan	7
GE Healthcare	8
Mediana	12
Medtronics	13
Midmark	14
Mindray	15
Nihon Kohden	16
Philips	17
Spacelabs	18
Zoe	19




Introduction

Welcome to the Exergen OEM Partner Vital Signs Monitor Catalog. The purpose of this catalog is to give Exergen Sales Reps and Distributors a clear understanding of which monitors are offered by all of the Exergen OEM Partners and which Exergen thermometers are connected to them, along with the connection type and any options for the cable replacements.

There are 2 types of cable connections, hard wired (HW) and quick release (QR). Only the quick release versions can have cables replaced. The hard wired versions have the cable permanently attached to the thermometer.

The list of the Quick Release (QR) cables offered by Manufacturer are below:

PART #	DESCRIPTION	MANUFACTURER	MONITOR	CABLE TYPE
124007	FRU CABLE 8' COILED – GE P/N 2044860-005	GE	Carescape V100	TTL (old TAT-5000-TTL version)
124075	FRU CABLE 8' COILED – GE P/N 2044860-018	GE	Carescape V100	TTL (new TAT-5000S-TTL version)
124208	CABLE QR REPLACEMENT	CAPSULE	SmartLinx Vitals Plus	RS-232
124299	CABLE QR REPLACEMENT	ZOE	740 Select	RS-232
124330	CABLE QR REPLACEMENT	DRAEGER	VISTA 120 SC	RS-232
		EDAN	IM-3	
		NIHON KOHDEN	SVM-7200	
		MEDIANA	V20	
		DIRECT SUPPLY	Attendant VSM	
		MEDTRONICS	RESPARRAY	
		ALFAMED	VITA i20	
124284	CABLE QR REPLACEMENT, MINDRAY	MINDRAY	Accutorr 7, VS-9, VS-900	RS-232




1	84001	ADHESIVE PRIMER CYANOACRYL
2	84001	ADHESIVE PRIMER CYANOACRYL

124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

Alfamed VITA i20

124234-AC-GN Deg C, Generic; Quick Release, RS-232 (CN)



124225-AC-QR Deg C, Generic; Quick Release, RS-232 (USA)



124208 – Cable, QR Replacement, CAPSULE

SmartLinux Vitals Plus

124239-AF-CT - TAT-5000S-RS232, ART, DEG F, Capsule Tech RS232, 4th Edition, quick release, RS-232




124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

Attendant® Touchscreen VSM

124234-AF-GN TAT-5000S-RS232-QR, ART – Deg F, Generic

124234-AC-GN Deg C, Generic; Quick Release, RS-232

124234-AC-GN Deg C, Gen




124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

IM3- Alfamed VITA i20

124234-AC-GN Deg C, Generic; Quick Release, RS-232 (CN)

124225-AC-QR Deg C, Generic; Quick Release, RS-232 (USA)




VSM Model: Portrait

Thermometer Model:

Exergen Item #:

TAT-5000S-USB-GE2-ART, DEG F	124227-AF-GE
TAT-5000S-USB-GE2-ART, DEG C	124227-AC-GE
TAT-5000S-USB-GE2-ORAL, DEG F	124227-OF-GE
TAT-5000S-USB-GE2-ORAL, DEG C	124227-OC-GE

USB connector, hard wired (HW)



VSM Model: GE Carescape V100

Thermometer Model:

Exergen Item #:

TAT-5000S-TTL,ART,DEG F,4th ed. GE P/N: 2044860-001, REV D	124003
TAT-5000S-TTL,ART,DEG C,4TH ED. GE P/N: 2044860-002, REV E	124004
TAT-5000S-TTL,ORAL,DEG C,4TH ED. GE P/N: 2044860-003, REV E	124005
TAT-5000S-TTL,ORAL,DEG F,4th ed GE P/N: 2044860-013, REV D	124006

15 pin connector, TTL quick release (QR)

124007 -FRU CABLE 8' COILED - GE P/N 2044860-005 (OLD TAT-5000-TTL)

124075 -FRU CABLE 8' COILED - GE P/N 2044860-018 (NEW TAT-5000S-TTL)



VSM Model: Corometrics 250C

Thermometer Model:

Exergen Item #:

TAT-5000S-CORO,ART,Deg F,4th ed GE P/N: 2036641-001, REV: F	124001
TAT-5000S-CORO,ART,Deg C,4th ed GE P/N: 2036641-002, REV: F	124002
TAT-5000S-CORO,ORAL,DEG F,4 th ed GE P/N: 2036641-003, REV: E	124010
TAT-5000S-CORO,ORAL,DEG C,4 th ed GE P/N: 2036641-004, REV: D	124011

15 pin connector, hard wired (HW), RS-232




VSM Model: GE CARESCAPE VC150 Vital Signs Monitor

Thermometer Model:

Exergen Item #:

TAT-5000S-USB,ART,DEG F,4th ed. GE P/N: 2044860-013, REV E	124260
TAT-5000S-USB,ART,DEG C,4th ed. GE P/N: 2044860-014, REV E	124261
TAT-5000S-USB-ORAL, DEG F, 4 TH ed. GE P/N: 2044860-015, REV E	124262
TAT-5000S-USB-ORAL, DEG C, 4 th ed. GE P/N: 2044860-016, REV E	124263

USB, hard wired (HW)




124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

V20

124234-AF-GN TAT-5000S-RS232-QR, ART – Deg F, Generic


124234-AC-GN Deg C, Generic; Quick Release, RS-232




124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

RespArray

124234-AC-GN Deg C, Generic; Quick Release, RS-232



124264-AF-MD - TAT-5000S-RS232-QR, ART DEG F, 3-200-0047, 4th edition



9	840081	ADHESIVE, PRI
8	819072	
7	840000	
6	819115	H
5	822120	HIRC
4	720830	SPOOL
3	840046	
2	822127	AUDI
1	819207	CO
ITEM NO. PART NO.		
UNLESS OTHERWISE SPECIFIED: DRAWN: J. LEE 291		
DIMENSIONS ARE IN INCHES CHECKED: G. TSOTSI 0		
TOLERANCES: APPD: J. LEE 02		
ANGULAR MATCH 1*		



124284 – CABLE, QR REPLACEMENT, MINDRAY – Ref: Mindray 040-003354-00

VS-9, VS-9000, Accutorr 7 – Quick Release, RS232, Hirose Plug

TAT-5000S-RS232-QR, MINDRAY	124237-AF-MR
TAT-5000S-RS232-MR,ART, DEG F, MINDRAY, CHINA	124233-AF-MR
TAT-5000S-RS232-MR,ART,DEG C, MINDRAY, CHINA	124233-AC-MR
TAT-5000S-RS232-MR, ORAL, DEG F, MINDRAY, CHINA	124233-OF-MR
TAT-5000S-RS232-MR, ORAL, DEG C, MINDRAY, CHINA	124233-OC-MR



SVM-7200



124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

SVM-7200

124234-AF-GN TAT-5000S-RS232-QR, ART – Deg F, Generic

124234-AC-GN Deg C, Generic; Quick Release, RS-232

PHILIPS


EarlyVue VS30



IntelliVue MP5SC



SureSigns VS4




Philips EarlyVue VS30, Philips IntelliVue MP5SC, Philips SureSigns VS4, RS232, Hard Wired, Hirose HR30-6 position plug

T5KS-RS232-PH, ART F, 4th ed. Philips P/N: 989803192931 C	124265
TAT-5000S-232EP, ART DEG F Philips P/N: 453564961221 A	124229-AF-EP



Qube®, Xprezzon®, Qube® Mini

124231-AF-GN - TAT-5000S-USB, ART, DEG F, 4TH ed, SPACELABS 010-2157-00, REV F - HARDWIRED



124299 CABLE, QR REPLACEMENT, REF ZOE MEDICAL P/N 421-0601

740 SELECT

124264-AF-MD - TAT-5000S-RS232-QR, ART, DEG F – MIDMARK PN:3-200-0047, 4TH ed , Quick Release