

EXERGEN OEM PARTNER VITAL SIGNS MONITOR CATALOG

Complete list of
VSM,
Thermometer,
and Cable
options with part
#'s

Table of Contents **Page #**

Introduction	2
GE Healthcare	3
Draeger	6
Nihon Kohden	7
Mediana	8
Capsule	9
Edan	10
Midmark	11
Mindray	12
Philips	14
Zoe	15
Spacelabs	16

Introduction

Welcome to the Exergen OEM Partner Vital Signs Monitor Catalog. The purpose of this catalog is to give Exergen Sales Reps and Distributors a clear understanding of which monitors are offered by all of the Exergen OEM Partners and which Exergen thermometers are connected to them, along with the connection type and any options for the cable replacements.

There are 2 types of cable connections, hard wired (HW) and quick release (QR). Only the quick release versions can have cables replaced. The hard wired versions have the cable permanently attached to the thermometer.

The list of the Quick Release (QR) cables offered by Manufacturer are below:

PART #	DESCRIPTION	MANUFACTURER	MONITOR	CABLE TYPE
124007	FRU CABLE 8' COILED – GE P/N 2044860-005	GE	Carescape V100	TTL (old TAT-5000-TTL version)
124075	FRU CABLE 8' COILED – GE P/N 2044860-018	GE	Carescape V100	TTL (new TAT-5000S-TTL version)
124208	CABLE QR REPLACEMENT	CAPSULE	SmartLinx Vitals Plus	RS-232
124299	CABLE QR REPLACEMENT	ZOE	740 Select	RS-232
124330	CABLE QR REPLACEMENT	DRAEGER	VISTA 120 SC	RS-232
124330	CABLE QR REPLACEMENT	NIHON KOHDEN	SVM-7200	RS-232
124330	CABLE QR REPLACEMENT	MEDIANA	V20	RS-232
124330	CABLE QR REPLACEMENT	EDAN	IM-3	RS-232
124284	CABLE QR REPLACEMENT, MINDRAY	MINDRAY	Accutorr 7, VS-9, VS-900	RS-232




VSM Model: Corometrics 250C

Thermometer Model:

Exergen Item #:

TAT5000S-CORO,ART, Deg F, 4th ed GE P/N: 2036641-001, REV: F	124001
TAT5000S-CORO,ART, Deg C, 4th ed GE P/N: 2036641-002, REV: F	124002
TAT5000S-CORO,ORAL, DEG F, 4 th ed GE P/N: 2036641-003, REV: E	124010
TAT5000S-CORO,ORAL, DEG C, 4 th ed GE P/N: 2036641-004, REV: D	124011

15 pin connector, hard wired (HW), RS-232



VSM Model: GE Carescape V100

Thermometer Model:

Exergen Item #:

TAT5000S-TTL,ART,DEG F,4th ed. GE P/N: 2044860-001, REV D	124003
TAT5000S-TTL,ART,DEG C,4TH ED. GE P/N: 2044860-002, REV E	124004
TAT5000S-TTL,ORAL,DEG C,4TH ED. GE P/N: 2044860-003, REV E	124005
TAT5000S-TTL,ORAL,DEG F,4th ed GE P/N: 2044860-013, REV D	124006

15 pin connector, TTL quick release (QR)

124007 -FRU CABLE 8' COILED - GE P/N 2044860-005 (OLD TAT-5000-TTL)

124075 -FRU CABLE 8' COILED - GE P/N 2044860-018 (NEW TAT-5000S-TTL)



VSM Model: GE CARESCAPE VC150 Vital Signs Monitor



Thermometer Model:

Exergen Item #:

TAT5000S-USB,ART,DEG F,4th ed. GE P/N: 2044860-013, REV E	124260
TAT-5000S-USB,ART,DEG C,4th ed GE P/N: 2044860-014, REV E	124261
TAT-5000S-USB-ORAL, DEG F, 4 TH ed. GE P/N: 2044860-015, REV E	124262
TAT-5000S-USB-ORAL, DEG C, 4 th ed. GE P/N: 2044860-016, REV E	124263

USB, hard wired (HW)

Dräger





124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

**VISTA 120 SC - 124234-AF-GN TAT5000S-RS232-QR, ART – Deg F, Generic;
124234-AC-GN Deg C, Generic; Quick Release, RS-232**




SVM-7200




124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

**SVM-7200 - 124234-AF-GN TAT5000S-RS232-QR, ART – Deg F, Generic;
 124234-AC-GN Deg C, Generic; Quick Release, RS-232**




124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

V20 - 124234-AF-GN TAT5000S-RS232-QR, ART – Deg F, Generic; 124234-AC-GN Deg C, Generic; Quick Release, RS-232



124208 – Cable, QR Replacement, CAPSULE

SmartLinx Vitals Plus - 124239-AF-CT - TAT5000S-RS232, ART, DEG F, Capsule Tech RS232, 4th Edition, quick release, RS-232




124330 CABLE, QR REPLACEMENT, TAT-5000S-RS232-QR

IM3-

124234-AC-GN Deg C, Generic; Quick Release, RS-232 (CN)



124225-AC-QR Deg C, Generic; Quick Release, RS-232 (USA)




124299 CABLE, QR REPLACEMENT, REF ZOE MEDICAL P/N 421-0601

Digital Vital Signs Device - 124264-AF-MD - TAT-5000S-RS232-QR, ART DEG F, 3-200-0047, 4th edition

mindray



SECTION A-A
SCALE 2:1



HIROSE PLUG PINOUT
VIEW FROM SOLDER SIDE

9	840081	ADHESIVE, PRI
8	819072	
7	840000	
6	819115	H
5	822120	HIRC
4	720830	SPOOL
3	840046	
2	822127	AUDI
1	819207	CO

ITEM NO. PARTNO
UNLESS OTHERWISE SPECIFIED: DRAWN: J. LEE SF1
DIMENSIONS ARE IN INCHES CHECKED: G. TIGOTI 0
TOLERANCES: ANGLE/ARE MACH 1*




124284 – CABLE, QR REPLACEMENT, MINDRAY – Ref: Mindray 040-003354-00

VS-9, VS-9000, Accutorr 7 – Quick Release, RS232, Hirose Plug

TAT-5000S-RS232-QR, MINDRAY	124237-AF-MR
TAT5000S-RS232-MR,ART, DEG F, MINDRAY, CHINA	124233-AF-MR
TAT5000S-RS232-MR,ART,DEG C, MINDRAY, CHINA	124233-AC-MR
TAT-5000S-RS232-MR, ORAL, DEG F, MINDRAY, CHINA	124233-OF-MR
TAT-5000S-RS232-MR, ORAL, DEG C, MINDRAY, CHINA	124233-OC-MR
CABLE, QR REPLACEMENT, MINDRAY	124284

PHILIPS




Philips EarlyVue VS30, Philips IntelliVue MP5SC, Philips SureSigns VS4, RS232, Hard Wired, Hirose HR30-6 position plug

T5KS-RS232-PH, ART F, 4th ed. Philips P/N: 989803192931 C	124265
TAT-5000S-232EP, ART DEG F Philips P/N: 453564961221 A	124229-AF-EP

ZOE


Medical



124299 CABLE, QR REPLACEMENT, REF ZOE MEDICAL P/N 421-0601

740 SELECT

124264-AF-MD - TAT-5000S-RS232-QR, ART, DEG F – MIDMARK PN:3-200-0047, 4TH ed , Quick Release



Qube®, Xprezzon®, Qube® Mini

**124231-AF-GN - TAT-5000S-USB, ART, DEG F, 4TH ed, SPACELABS 010-2157-00, REV F -
HARDWIRED**